

Mammalian Cell Transfection

Date _____

Cell line for the transfection _____

- Inoculate a 100mm culture dish with exponentially growing cells at 5×10^5 cells per dish no more than 24 hours before the transfection (for L293 cells, pass 1:6 at the end of the day for transfection the next day).
- The transfection should be done at the end of the day so that they will sit in the CaCl_2 and BBS only overnight. Follow the steps below to set up the transfection.
 - Dilute the desired amount of DNA with distilled, deionized, sterile water to 450 μL . (for L293 cells, 10 μg DNA works well). **USE ONLY POLYSTYRENE TUBES !!!**
 - Add 50 μL of the sterile 2.5 M CaCl_2 solution.
 - Set a Vortex mixer to #3 and while vortexing the tube, slowly add 500 μL sterile 2XBBS.
 - Allow the mixture to incubate for 20 minutes at room temperature.
 - Gently mix the precipitate to ensure adequate suspension of the precipitate. Then add the suspension dropwise to the media while swirling the plate to ensure even distribution.
 - After 12-24 hours (overnight), gently remove the media and wash the cells twice with PBS (without Ca or Mg) and add fresh complete media.
- For transient transfections, the cells may be collected 48 hours after the transfection.

2.5 M CaCl_2

18.37 g

Bring up to 50 mL w/ddH₂O

filter sterilize

2 X BBS

50 mM BES

280 mM NaCl

1.5 mM Na_2HPO_4

pH to exactly 6.95 with HCl

filter sterilize

for 500 mL

5.33 g

8.18 g

0.1065 g

Tube #	DNA Description	Amount of DNA for _____ μg	Water	2.5 M CaCl_2	2 X BBS
1				50 μL	500 μL
2				50 μL	500 μL
3				50 μL	500 μL
4				50 μL	500 μL
5				50 μL	500 μL
6				50 μL	500 μL
7				50 μL	500 μL
8				50 μL	500 μL
9				50 μL	500 μL
10				50 μL	500 μL
11				50 μL	500 μL
12				50 μL	500 μL
13				50 μL	500 μL
14				50 μL	500 μL
15				50 μL	500 μL
16				50 μL	500 μL
17				50 μL	500 μL
18				50 μL	500 μL
19				50 μL	500 μL
20				50 μL	500 μL