

Silver Staining of proteins in a Gel

Stock Solutions

50 % Methanol; 10 % Acetic Acid

10 % Methanol; 5 % Acetic Acid

100 mM DTT

Farmer's Reducer

Solution A - Add 37.5 g potassium ferricyanide and bring up to 500 ml with dH₂O.

Solution B - Add 120 g anhydrous sodium thiosulfate and bring up to 500 ml with dH₂O.

To be prepared fresh (Prepare only what you will use)

0.012 M AgNO₃ -- Add 0.041 g AgNO₃ to 20 ml dH₂O.

0.28 M Na₂CO₃ + Formaldehyde -- Add 0.6 g Na₂CO₃ and 10µl HCHO to 20 ml dH₂O.

Procedure

1. Fix the gel in 50 % methanol: 10 % Acetic acid for 15 minutes (***30 minutes***)* or overnight.
2. Wash in 10 % methanol: 5 % Acetic acid for 15 minutes (***30 minutes***)*.
3. Wash 1 minute (***5 minutes***)* in dH₂O. Repeat 3-5 times
4. Add 2 µl of 100 mM DTT to 100 ml dH₂O add soak for 15 minutes (***30 minutes***)*.
5. Incubate in 0.012 M AgNO₃ for 15 minutes (***30 minutes***)*.
6. Wash twice (5 seconds each) in dH₂O.
7. Rinse briefly in approx. 25 ml Na₂CO₃/HCHO.
8. Repeat rinse.
9. Add remaining Na₂CO₃/HCHO and rock gently until bands begin to appear.
10. Remove Na₂CO₃/HCHO when the gel is nearing the desired intensity, and wash repeatedly in large volumes of dH₂O.

Discard all washes & solutions in aqueous waste container NOT IN SINK!!

To Reduce Background Staining

1. Mix 1 ml Farmers Solution A with 4 ml Farmers Solution B and bring up to 100ml with dH₂O.
2. Add the Farmers solution to the gel and watch closely until the gel is nearing the desired destaining intensity.
3. Wash repeatedly in large volumes of dH₂O.
4. Some destaining may occur during the dH₂O washes, so it is best to start the washing a little earlier.

Discard all washes & solutions in aqueous waste container NOT IN SINK!!

* Incubation times appearing in parenthesis, bold, and italic are to be used for higher percentage acrylamide gels.

Silver Staining of proteins on a nitrocellulose blot

OBVIOUSLY NOT FOR USE WITH A BLOT THAT HAS BEEN BLOCKED!!

1. Wash blot in sterile distilled water and leave in water until ready.
2. Make the following solutions:
 - a. 40 % Na Citrate
 - b. 20 % Silver nitrate
 - c. 20 % $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ (to be made up fresh every time!!)
3. To a 15 ml tube, add 9 ml water
4. Add 0.5 ml 40 % Na Citrate and mix.
5. Add 0.4 ml 20 % $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ and mix.
6. Add 0.1 ml 20 % Silver nitrate and mix very well!! It will form a dark precipitate but this will go back into solution quickly. Use immediately !!!
7. Remove water from blot and quickly add staining solution.
8. Mix vigorously to avoid over staining only parts of the blot.
9. Staining is very rapid so as soon as the bands show up, remove the stain and immediately wash with plenty of running water (place blot directly under stream of water).